

## Vista Analytical Laboratory CARB Method 428

CDD/CDF	QL (pg/sample)
2,3,7,8-TCDD	10.0
2,3,7,8-TCDF	10.0
1,2,3,7,8-PeCDD	50.0
1,2,3,7,8-PeCDF	50.0
2,3,4,7,8-PeCDF	50.0
1,2,3,4,7,8-HxCDD	50.0
1,2,3,6,7,8-HxCDD	50.0
1,2,3,7,8,9-HxCDD	50.0
1,2,3,4,7,8-HxCDF	50.0
1,2,3,6,7,8-HxCDF	50.0
1,2,3,7,8,9-HxCDF	50.0
2,3,4,6,7,8-HxCDF	50.0
1,2,3,4,6,7,8-HpCDD	50.0
1,2,3,4,6,7,8-HpCDF	50.0
1,2,3,4,7,8,9-HpCDF	50.0
OCDD	100.0
OCDF	100.0