

LCMS Bottle Order Form

Vista Analytical Laboratory - PFAS Bottle Order Form

Client/Company Name:	<input type="text"/>	Date to be Received:	<input type="text"/>
Ship To Address:	<input type="text"/>	*Time to be Received:	<input type="text"/>
		**Client's FedEx account:	<input type="text"/>
Project Name:	<input type="text"/>	Requester Name:	<input type="text"/>
Project Number:	<input type="text"/>	Requester Email:	<input type="text"/>
Recipient Name:	<input type="text"/>	Requester Phone:	<input type="text"/>
Recipient Email:	<input type="text"/>	Recipient Phone:	<input type="text"/>

* Note if required by 8:30 am or 10:30 am, otherwise bottles are shipped via ground.
 ** FedEx Account Number required for rush samples, and for 8:30 am or 10:30 am deliveries.

Unmodified EPA 537.1 (Drinking Water) / PFAS by Isotope Dilution Bottleware

Sample Container Type	Number of Containers	
125 mL HDPE bottle	<input type="text"/>	250 mL HDPE bottle <input type="text"/>
125 mL HDPE with Trizma	<input type="text"/>	250 mL Polypropylene w/Trizma <input type="text"/>
6 oz. HDPE jar (solid)	<input type="text"/>	250 mL HDPE with Trizma

2 bottles per AQ sample (6 bottles per samples designated MS/MSD)
 1 jar per soil/solid sample (including MS/MSD)

Ultra Pure Water

Ultra Pure Water	Number of Bottles		Cooler Size	Number of Coolers	
	w/water	empty		w/ water	empty
125 mL HDPE	<input type="text"/>	<input type="text"/>	Full Size	<input type="text"/>	
125 mL HDPE with Trizma	<input type="text"/>	<input type="text"/>	Small	<input type="text"/>	
250 mL HDPE	<input type="text"/>	<input type="text"/>	250 mL HDPE with Trizma		
250 mL Polypropylene w/Trizma	<input type="text"/>	<input type="text"/>			
1 Liter HDPE	<input type="text"/>				
1 gallon HDPE	<input type="text"/>				

Each FRB requires: 2 bottles of water and 2 empty bottles.

Standard Supplies: Sample Container Labels, CoCs, Temp Blanks, Custody Seals.

Client Notes: